

Leica RM2135 rotary microtome

The versatile easy-to-use rotary microtome



Functionality, safety and comfort

Future-based concepts must offer solutions to the multiple challenges of practical work. Leica's rotary microtomes - based on 125 years of experience - meet these demands perfectly. The new Leica RM2135 features advanced technology, extraordinary safety standards along with an ergonomic design. The instrument presents technical innovations and improvements such as a patented gravitational force compensation system; a coarse feed wheel located in a convenient position to ensure stress-free working; cross roller bearings for the vertical and horizontal specimen feed to guarantee high-quality sections - features designed to meet even the highest standards in paraffin and hard-sectioning techniques.




Precision: Specimen orientation *

An easy-to-use specimen orientation system with "anti-tilt feature" ensures accurate and simple orientation of the specimen. A vertical stroke length of 70 mm guarantees a maximum sectioning range of 69 mm. This enables sectioning of large specimens.



Efficiency: The new quick release specimen clamping system *

The new guick release clamping system can be operated with one hand and enables extremely fast changing of specimen clamps.



The new section waste tray is integrated in the microtome base plate, offering a great advantage

to the user when removing section debris. Sectioning waste is directly transferred from the knife into the waste tray. The volume can be expanded when using disposable waste bags.





Cross roller bearings for utmost precision:

The cross roller bearings for vertical and horizontal specimen feed ensure accurate reproducibility of the section thickness. Even when using high cutting speeds and strong cutting forces the stroke will be free from vibrations. Leica is the first manufacturer to use this technology in manually operated rotary



Unsurpassed flexibility: The new knife holder system *

The universal knife holder base with lateral displacement function is suitable for holding different knife holders. The knife holder can be moved laterally to take full advantage of the whole cutting length without having to readjust the clearance angle.




The patented gravitational force compensation system

Innovation:

By means of a pre-stressed spring, the patented lever system offsets the centrifugal forces that arise during the sectioning process. Thus, there is no need for mounting an approximately 2 kg counterweight inside the handwheel. The new system ensures a very smooth and easy turn of the handwheel to obtain a uniform stroke for excellent, perfectly stretched sections.

. Standard: Counterweight inside the handwheel



Inventive: The new section waste tray *

Stop: Safety locking

The quick locking lever locks the handwheel quickly and safely in any position. The new 12 o'clock locking system automatically locks the handwheel in the upper position by simply pulling the handle outwards.

Personalized comfort: The coarse feed wheel

The coarse feed wheel is ergonomically positioned to ensure stress-free work conditions. You have the option to order the unit with the coarse feed wheel rotating clockwise or counterclockwise. The (optional) mechanical trimming device is placed very close to the coarse feed wheel. As an alternative, the RM2135 can also be ordered with the coarse feed wheel located at the front of the microtome base plate.







Leica RM2135 – Instrument features

- Ergonomically designed
- Mechanical high-precision specimen feed
- Backlash-free micrometer feed system in dust-proof housing
- Patented mass-balanced system to compensate gravitational forces
- · Smoothly running handwheel with quick locking system
- Specimen retraction can be disengaged
- Coarse feed wheel in front or on the side
- Totally enclosed housing
- Easy to service

- Low-maintenance microtome with cross roller bearings for vertical and horizontal specimen movement
- Quick release clamping system *
- Simple changing of specimens
- Mechanical trimming function; at two steps (10 μm and 30 $\mu m)^*$
- Specimen orientation with anti-tilt feature *
- Ergonomic hand rest *
- Section waste tray with holder for disposable waste bags *
- Knife holder base with lateral displacement feature *

* optional accessories

Leica RM2135 – Technical specifications

Section thickness setting:	1 to 60 μm
	from 1 to 10 μm in steps of 1 μm
	from 10 to 20 μm in steps of 2 μm
	from 20 to 60 μm in steps of 5 μm
Section thickness indication:	display window
Coarse advance:	manually via coarse feed wheel
Horizontal specimen feed:	
Vertical specimen stroke lengt	h:
Sectioning range without retra	ction:69 mm (at 1 μm)
Sectioning range with retraction	on:62 mm


Specimen orientation:	
Dimensions (L x W x H):	525 x 310 x 295 mm
Weight (without accessories):	

Technical specifications subject to change. Wide range of accessories on request.

Up-to-date development, production and quality control procedures certified under DIN EN ISO 9001 ensure highest quality and reliability.

For further information please contact your local Leica sales company

Australia:	Gladesville/NSW	Phone +61 2 9897 9700	Fax +61 2 9897 8358
Austria:	Vienna	Phone +43 1 486 80 50	Fax +43 1 486 80 50 30
Canada:	Richmond Hill/Ontario	Phone +1 905 762 2000	Fax +1 905 762 8937
Denmark:	Herlev	Phone +45 4454 0101	Fax +45 4454 0111
France:	Rueil-Malmaison	Phone +33 1 4732 8585	Fax +33 1 4732 8586
Germany:	Bensheim	Phone +49 6251 1360	Fax +49 6251 136 155
Italy:	Milan	Phone +39 02 57 486.1	Fax +39 02 57 40 3273
Japan:	Tokyo	Phone +81 3 5435 9603	Fax +81 3 5435 9615
Korea:	Seoul	Phone +82 2 514 6543	Fax +82 2 514 6548
Netherlands:	Rijswijk	Phone +31 70 4132130	Fax +31 70 4132139
Portugal:	Lisbon	Phone +351 1 388 9112	Fax +351 1 3854 668
Republic of China:	Hong Kong	Phone +852 2564 6699	Fax +852 2564 4163
Singapore:	Singapore	Phone +65 67 79 78 23	Fax +65 67 73 06 28
Spain:	Barcelona	Phone +34 93 494 9530	Fax +34 93 494 9532
Sweden:	Sollentuna	Phone +46 8 625 4545	Fax +46 8 625 4510
Switzerland:	Glattbrugg	Phone +41 1 809 34 34	Fax +41 1 809 34 44
United Kingdom:	Milton Keynes	Phone +44 1 908 246 246	Fax +44 1 908 609 992
USA:	Bannockburn/Illinois	Phone +1 847 405 0123	Fax +1 847 405 0164



...and more than 100 national distributors.

Leica Microsystems Nussloch GmbH Heidelberger Strasse 17-19 D-69226 Nussloch Telephone: (06224) 143-0 Fax: (06224) 143 200 e-mail: histo_info@leica-microsystems.com www.leica-microsystems.com

