

# Comparison Chart

## Dakewe CT520 + Leica CM1950

**RANKIN**

1-877-882-3679



**DAKEWE**

**Leica**  
MICROSYSTEMS



Cools to -35C (Chamber)	✓	✓
Specimen Head Cools to -50C	✓	✓
Left-hand control buttons	✓	✓
Built-in heat extractor	✓	✓
Built-In anti roll plate	✓	✓
Accepts Leica cryo chucks	✓	✓
Peltier on freezing shelf	✓	✓
2-in-1 Blade Holder	✓	✓
UV Disinfection	✓	✓
Numbered freezing shelf	✓	✓
Blade Holder Temp Control	✓	✗
Intuitive Color Touchscreen	✓	✗
Slide placement area	✓	✗
Scheduled UV cleaning	✓	✗
Drying mode	✓	✗
On-screen visual coarse feed location	✓	✗
2-Year Warranty	✓	✗
Motorized Sectioning and/or Vacuum	✗	✓ *optional

