

# GL119



The pathology work station provides side downward exhaust to the dissecting area. The unit can be operated from either side.

## STANDARD DESIGN FEATURES

Length: .....Per Customer Requirements  
Height/Overall: .....40 in./101 cm.  
Width: .....26 in./66 cm.

- All stainless steel construction of the entire work area and all exterior closure panels.
- Exhaust duct for connection to building system
- Downward exhaust on both sides of dissecting area
- Integral sink with mixing faucet
- Hot and cold hand spray assembly

- Heavy-duty commercial waste disposer, 1/2 hp and switch
- (4) Instrument drawers
- Pull out writing platform
- Polyethylene dissecting board (Photo Blue)
- GFCI duplex receptacle w/waterproof cover

All measurements in inches/centimeters

We also carry Biopsy Bags



Call for more information  
(800) 410-0100

## MODEL GL119 SIDE EXHAUST PATHOLOGY WORKSTATION

## Optional Accessories

- Heavy-duty commercial disposer (3/4 hp, 1 hp, 2 hp)
- Eyewash drench assembly
- Eyewash assembly, faucet mounted
- CO2 valves/vacuum valves/ compressed air valves
- Hands free soap dispensers
- Specimen rinse
- Digital platform scale, 4 lbs./2 kg capacity, economy or premium grades available
- Removable measuring rule
- Computer monitor and keyboard stand premium or standard grades
- CPU wall mount
- Digital voice recorder
- Articulating arm fluorescent light w/magnifier
- Articulating arm examination light w/magnifier
- Flex arm incandescent light w/magnifier
- Flex arm halogen lights
- Seismic Anchors

### Options based on available space

When planning an order, we will assist you in selecting options to complete your system. Please feel free to contact us.

Custom layout and designs are available.

