

Preliminary Specifications

# Signature<sup>®</sup> Cassette Printer

Print directly onto cassettes for more efficient and more accurate processing of lab specimens.

Primera's Signature Cassette Printer can significantly increase the efficiency of your lab while helping to reduce the risk of misidentification of specimens. You can print text, graphics and logos along with linear and 2D bar codes directly onto the face of cassettes, making your workflow more efficient while increasing patient safety.

## Key Features and Benefits:

- » On-demand or batch mode printing
- » Black or color printing
- » 4 cassette hoppers
- » 40 cassettes per hopper
- » Up to 11 cassettes per minute
- » Thermal transfer printing at 300 dpi
- » Ribbon types: Black – 5,000 prints / CMYK – 1,000 prints
- » Prints on standard cassettes with 35° or 45° angles with either lid on or lid off, no adjustments needed
- » Output Capacity: 7 cassettes / 17 with extension tray
- » Chemical-resistant ink – ensures reliable identification of cassettes
- » USB interface – ability to integrate with Laboratory Information Systems
- » Available Q1 2014



**PRIMERA**  
TECHNOLOGY, INC.

Concept Drawing

510748-080613