



**SlideID**

**Complete  
Printer  
Solution**

# MICROSCOPE SLIDE LABEL PRINTING

the EASY way...



## READY-TO-PRINT

*SlideID* is a complete, ready-to-print microscope slide label printer package that comes with every thing you need to print histology grade slide labels—printer, software, labels and ribbon.

Take the guesswork out of it, and do it right the first time. SlideID is the easy, efficient and economical way to eliminate manual, hand-labeling of microscope slides and implement laboratory best practices.



## RESISTANT TO STAINS, SOLVENTS, AND XYLENE

*SlideID* microscope slide labels are applied BEFORE the staining process. The labels and thermal transfer ribbon are optimized for harsh chemicals and will not absorb eosin or hematoxylin coloring.



## SCALABLE

*SlideID* printer packages are ideal for the small to mid-sized lab just starting with slide printing, as well as high-volume labs that need multiple printers at each workstation. They function with, or without LIS. SlideID printers can be setup for batch printing or with on-demand printing for "laboratory best practices."

# SlideID

## Complete Printer Solution



### SlideID Complete Printer Solutions

#### BASIC PACKAGE:

SlideID CL-E331 Slide Label Printer by Citizen Systems-300dpi Thermal Transfer Label Printer setup for Microscope Slide Label Printing

**SlideID Slide Printing Software** - Simple, Microscope Slide Labeling Software with templates for pathology best practices and free text, keyboard data entry.

**SlideID Microscope Slide Labels** — 1,800 per roll. Histology grade. Applied BEFORE staining. Xylene and Solvent Resistant. Do not absorb hematoxylin or eosin coloring.

**SlideID Thermal Transfer Ribbon**— Special resin ribbon optimized for microscope slide label printing. Solvent and xylene resistant. 984' Prints approximately 12,000 +/- microscope slide labels.

#### PRO PACKAGE (includes all of above):

**Denso AT30Q Barcode Scanner** – Optimized for healthcare applications. Reads difficult and small 2D barcodes. Pre-programmed for your internal barcode specifications.

#### PRO WORKSTATION (includes all of above):

**All-in-One PC. or Tablet PC** preconfigured with SlideID software and Denso barcode scanner. Arrives ready-to-print setup to your internal specifications.

